

Technical Data Sheet


Comparing Lanthanum Hexaboride (LaB₆) and Cerium Hexaboride (CeB₆) Cathodes

#80920 - 80933

The performance and lifetime of the hexaboride cathode are determined by several factors: vacuum level, cathode temperature, impurity level, crystal orientation, tip shape, and mount design.

CeB₆ exhibits a lower evaporation rate


CeB₆ has an evaporation rate at normal operating temperatures near 1800 K that is lower than that of LaB₆. So long as care is taken to operate the cathode below 1850 K, CeB₆ should maintain an optimum tip shape longer, and therefore last longer.




Evaporation Rate Comparison for LaB₆ and CeBixTM at 2 x 10⁻⁸ Torr

CeB₆ is more resistant to carbon contamination

In laboratory tests, CeB₆ has proven to be more resistant to the negative impact of carbon contamination than LaB₆. This increases stability and results in less time to reach stable operation and more rapid recovery from contamination events.




Recovery from Carbon Contamination Comparison between LaB₆ and CeBixTM




Carbon Contamination Comparison between LaB₆ and CeBix[™]

Reduced effects of O₂ on CeBix[™]

The effects on oxygen on current density is shown in the figures below for LaB₆ and CeBix[™].



CeBix™



Material Data for LaB₆ and CeB₆

Parameter	Units	LaB ₆	CeB ₆
Stoichiometry	N/A		~6
Metal Impurities	ppm by wt.		<30
Density	g/cm ³	4.72	4.80
Coefficient of Thermal Expansion	$\alpha \times 10^6$	5.6	6.2
Electrical Resistivity	$\mu\Omega\text{-cm}$	~50	~65
Effective Work Function (100) at 1800 K	eV	2.70	2.65
Spectral Emissivity at 0.65 microns	N/A	0.765	0.779
Evaporation rate at 1800 K (UHV)	g/cm ² /s	2.2×10^{-9}	1.6×10^{-9}
Orientation limit for specific orientation	degrees		<2
Pyrolytic block mount resistance @ 1800 K	Ohms		1.45